

BASIC EMS RESCUE BACKPACK

Basic "14030" EMS Rescue Backpack

Description: Incorporates advanced features of modern outdoor long-haul hiking packs with the realistic needs and governmental requirements of Military SWAT, and special Ops medics. Created to help Tactical Medics carry necessary medical equipment in special operation situations.

- Foam-lined construction maintains pack shape and protects contents
- Main compartment houses three removable sub-packs that are held in place with hook-&-loop backing
- Internal sub-packs feature durable urethane windows for quick identification and access to equipment
- Quick-zip access to main compartment with reverse coil zippers to offer fast, smooth action while providing protection from dirt and wear
- Single-buckle top closure allows quick in and out
- Ergonomic, padded shoulder straps, back panel and waistbelt offer comfort during long transport
- padded, breathable back panel holds hydration bladder
- Large side pockets allow quick access to vital sign tools
- Whistle built into sternum-strap buckle adds safety
- Load handle for quick/easy carry
- Large zipper pulls allow easy, quick access even while wearing gloves
- Dimensions: 15" x 21" x 9.5"
- Packed 2/case

