

Folding Stretcher, Aluminum

www.basic-medicalsupply.com

Folding Stretcher, Aluminum

Designed for an easy stowage, this stretcher can be easily overlapped on other stretchers of the same kind to store a high number of them in a small space. Light and easy to handle, not foldable-Frame in sturdy aluminum. Complete with 2 quick unhooking belts for the immobilization of the patient-Orange PVC sheet fire proof, water proof, anti-tear and easy to clean.

Product Size(L*W*H):1860*480*180MM

Carton Size:1880*500*200mm

Load Limit:170kg

N.W.:5kg G.W.:11kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-FL001	Orange, Size: 1860*480*180mm, Material: Aluminum	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com