

Folding Stretcher, Aluminum

Folding Stretcher, Aluminum

Folding stretcher strengthened by two steel bars that allow its use even for heavy transports. It is foldaway for an easy location. Blue, high quality oxford sheet is characterized for its fire proof, and water proof, anti-tear and easy to clean.

Product Size (L*W*H):1850*520*18mm

Folded Size (L*W*H):920*520*90mm

Carton Size (2pcs): 960*530*160mm

Load Limit: 159kg

N.W.:6kg G.W.: 13kg

ITEM #	DESCRIPTION	QTY	PRICE
MS-FL005	Blue, Size: 1850*520*18mm, Material: aluminum	Each	

