

Basic Medical Supply (CHINA) 11-1703, # 225 Jinxiang Road Pudong, Shanghai, China Tel: +86 (21) 3175 0988 Fax: +86 (21) 3175 0989 **Basic Medical Supply (USA)** 19902 Flax Flower Dr. Richmond, TX 77407 Tel:+1-832-535-3107 Fax:++1-832-535-3107

Easy Fold Wheeled Stretcher,w/2 straps

Easy Fold Wheeled Stretcher,w/2 straps

Made of 1" sturdy anodized aluminum tubing frame with lush folding hinges. Four inch wheels and posts fold easily out of the way for compact storage. Cover is fabricated from tough 18 oz. green vinyl coated nylon. Cover is fungus and rot resistant and easily wiped clean. 2" wide patient restraint straps with automotive seatbelt type buckles are standard equipment.Dimensions: Open: 74" L x 20-1/2" W x 9" H Folded: 36-1/2" L x 20-1/2" W x 3" H

ITEM #	DESCRIPTION	QTY	PRICE
MS-FL006	Green, Size: 74"Lx20-1/2"Wx9"H, 2pcs/Case	Each	
MS-FL106	Orange, Size:74"Lx20-1/2"Wx9"H, 2pcs/Case	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com



Basic Medical Supply (CHINA) 11-1703, # 225 Jinxiang Road Pudong, Shanghai, China Tel: +86 (21) 3175 0988 Fax: +86 (21) 3175 0989 Basic Medical Supply (USA) 19902 Flax Flower Dr. Richmond, TX 77407 Tel:+1-832-535-3107 Fax:++1-832-535-3107