

Basic Medical Supply (CHINA) 11-1703, # 225 Jinxiang Road Pudong, Shanghai, China Tel: +86 (21) 3175 0988 Fax: +86 (21) 3175 0989

Basic Medical Supply (USA) 19902 Flax Flower Dr. Richmond, TX 77407 Tel:+1-832-535-3107 Fax:++1-832-535-3107

www.basic-medicalsupply.com

Portable transport unit

Portable transport unit

The portable transport unit is used to transport, transfer or rescue patients from areas inaccessible to stretchers.

It can then be used to transfer the patient from the EMS stretcher to the hospital bed after transport and folds down into a small package weighing only one pound.

The portable transport unit allows up to six rescuers to assist in a patient's transport or transfer by featuring 12 handles-two more than the original model that was introduced last year.

The new design also features an additional six inches in width and now measures a generous 40 by 80 inches. Its durability has been improved to stand up to the most rigorous conditions and ensure safe handling of patients in challenging environments.

Specifications: 12 handles, 6 per side fluid-resistant rugged non-woven polypropylene spun-bond construction easily foldable weighs less than one pound transports up to 1,500 pounds

ITEM #	DESCRIPTION	QTY	PRICE
22021	Grey; 200x100cm; Non-woven with PVC coat	Each	



info@basic-medicalsupply.com

www.basic-medicalsupply.com

